

**IMMUNOHISTOCHEMISTRY PROTOCOL: DLX-2 (TES-1)** 08.93

D.D. Eisenstat, MD; A. Bulfone, MD; M. H. Porteus, MD, PhD; J.L.R. Rubenstein, MD, PhD\_

**General:****1. Frozen embryo (-80 C) sectioned in cryostat (-20 C), 10 microns thickness**

Use Fisher Superfrost Plus slides at room temperature  
Once section is dry, place in slide box on dry ice  
When sectioning complete, slide boxes placed under vacuum, slightly ajar,  
at 4 C for at least 2 hours, then 3-4 desiccant capsules (Humicap 386)  
added, box placed in freezer bag with desiccant crystals (Drierite #23001  
anhydrous calcium sulfate) at -80 C.

**2. Stain reference slides (every 10 sections) with cresyl violet 0.5%**

(Chroma brand, filtered before use.)

- (a) Slides to room temperature, allow condensation to thaw
- (b) Fix in 4% paraformaldehyde (freshly prepared) in 1X PBS pH 7.5  
for 10 min at RT
- (c) Wash 3 times with water
- (d) Stain in cresyl violet 1-5 min at RT
- (e) Wash 3 times in water
- (f) Dehydrate in graded ethanols, 2 min each:  
50%, 75%, 85%, 95%, 100%, 100%, xylene, xylene
- (g) Mount with Permount (Fisher); place on coverslip first (to be done in  
hood)

**[20X PBS in 0.5 liter:** 10 g Na<sub>2</sub>HPO<sub>4</sub> (140mM), 3.6 g NaH<sub>2</sub>PO<sub>4</sub> (60 mM),  
76 g NaCl (2.6 M); pH to exactly 7.5; autoclave]

**[4% paraformaldehyde in 0.25 liter:** 10 g paraformaldehyde  
(Polysciences #00380), 12.5 mL 20x PBS pH 7.5, dH<sub>2</sub>O to 250 mL;  
dissolve until clear at 65 C, filter (0.2µM Nalgene) in hood, bring to room  
temperature.]

**IMMUNOHISTOCHEMISTRY PROTOCOL: DLX-2 (TES-1)**

08.93

**Protocol A: VECTOR *Elite* ABC Kit**

(modified by D. Eisenstat, MD)

- 1) Primary Antibody to Dlx-2 is **M650**, dilution 1:50 ( rabbit anti-mouse)  
[ 15  $\mu$ L blocking serum, 3  $\mu$ L antibody, 7.5  $\mu$ L 20X PBS pH 7.5, 7.5  $\mu$ L 1% thimerosal (Sigma #T5125), 117  $\mu$ L dH<sub>2</sub>O: total volume 150  $\mu$ L]
- 2) Secondary Antibody (**Vector BA-1000/ PK6101 kit**) biotinylated goat anti-rabbit, dilution 1:200  
[ 150  $\mu$ L stock blocking serum, 10 mL 1X PBS , 50  $\mu$ L biotinylated antibody stock]
- 3) Blocking serum (**Vector S-1000**): goat
- 4) Control antibody (**Vector I-1000**): rabbit IgG, dilution 1:1250 or 1:2500
- 5) Tertiary substance: (**Vector PK 6100/ 6101**) Elite ABC  
[2 drops (100  $\mu$ L) solution 'A' to 5 mL 1X PBS, mix, add 2 drops 'B', mix well]
- 6) Chromogenic substrate: DAB (**Vector SK-4100**)  
[ 2 drops buffer stock solution added to 5 mL water, mix; add 4 drops DAB stock solution, mix; add 2 drops hydrogen peroxide solution, mix.]

**Day 1**

1. Thaw frozen sections at RT
2. Fix in 4% paraformaldehyde at RT 15-20 min
3. Wash 3 times in 1X PBS pH 7.5, 5 min/wash at RT
4. Add water to humidifier chamber
5. Dry excess PBS off slides; do not let slides dry out
6. Circle tissue section with Lipshaw ISOLATOR pen
7. Add 150  $\mu$ L goat serum, undiluted.
8. Incubate overnight at 4 C.

**Day 2**

9. Prepare primary antiserum or control antibody in 10% blocking serum, 0.05% thimerosal ,1X PBS pH 7.5
10. Blot off excess serum and add 150 $\mu$ L primary antibody
11. Incubate 1-2 hours at RT
12. Wash as in step 3.

13. Prepare secondary antibody . Add 150  $\mu$ L per slide.
14. Incubate 30-60 minutes at RT
15. Blot excess. Wash as per step 3.

## **IMMUNOHISTOCHEMISTRY PROTOCOL: DLX-2 (TES-1)**

08.93

### **Protocol A: VECTOR *Elite* ABC Kit**

#### **Day 2 (continued)**

- 16.\* Prepare 250 mL of 3%# hydrogen peroxide (Sigma) in 100% methanol (Fisher). Quench endogenous peroxidases by washing for 30 min at RT. (\* optional step; use if high background level of endogenous peroxidases) (# try 0.3% H<sub>2</sub>O<sub>2</sub> first)
17. Blot excess. Wash as in step 3.
18. Prepare tertiary molecule: let stand at RT 30 min prior to use
19. Add 150  $\mu$ L to slide. Incubate 30 min at RT.
20. Blot excess. Wash as in 3.
21. Prepare DAB substrate. Waste and excess DAB treated with 2% bleach prior to discarding.
22. Add 150  $\mu$ L . Allow staining for 1- 30 minutes; usually 2-10 min. Observe under microscope.
23. Stop reaction in water for 5 min.
24. Dehydration in graded ethanols then xylene. (refer to cresyl violet protocol).
25. Mount with permount.

**IMMUNOHISTOCHEMISTRY PROTOCOL: DLX-2 (TES-1)**

08.93

**Protocol B: Biotin / Streptavidin / Peroxidase**

- 1) Primary Antibody to Dlx-2 is **M650** dilution 1:50 ( rabbit anti-mouse)  
[ 15  $\mu$ L blocking serum, 3  $\mu$ L antibody, 7.5  $\mu$ L 20X PBS pH 7.5, 7.5  $\mu$ L 1% thimerosal (Sigma #T5125), 117  $\mu$ L dH<sub>2</sub>O: total volume 150  $\mu$ L]
- 2) Secondary Antibody (**Boehringer Mannheim (BMB) #1214659**) Sheep anti-rabbit IgG biotinylated F(ab)' fragment, dilution 1:200  
[ 150  $\mu$ L blocking serum, 7.5  $\mu$ L antibody, 75  $\mu$ L 20X PBS pH 7.5, 1267.5 $\mu$ L dH<sub>2</sub>O: total volume 1500  $\mu$ L]
- 3) Blocking serum (**Sigma S-7773**): sheep
- 4) Control antibody (**Vector I-1000**): rabbit IgG, dilution 1:1250 or 1:2500
- 5) Tertiary substance: (**BMB #1089-153**) Streptavidin-peroxidase, dilution 1:1000  
[ 1.5  $\mu$ L streptavidin-peroxidase, 75  $\mu$ L 20X PBS pH 7.5 , 1423.5  $\mu$ L dH<sub>2</sub>O]
- 6) Chromogenic substrate: DAB (**Vector SK-4100**)  
[2 drops buffer stock solution added to 5 mL water, mix; add 4 drops DAB stock solution, mix; add 2 drops hydrogen peroxide solution, mix. ]

**Day 1**

1. Thaw frozen sections at RT
2. Fix in 4% paraformaldehyde at RT 15-20 min
3. Wash 3 times in 1X PBS pH 7.5, 5 min/wash at RT
4. Add water to humidifier chamber
5. Dry excess PBS off slides; do not let slides dry out
6. Circle tissue section with Lipshaw ISOLATOR pen
7. Add 150  $\mu$ L sheep serum, undiluted.
8. Incubate overnight at 4 C.

**Day 2**

9. Prepare primary antiserum or control antibody in 10% blocking serum, 0.05% thimerosal ,1X PBS pH 7.5
10. Blot off excess serum and add 150 $\mu$ L primary antibody
11. Incubate 1-2 hours at RT
12. Wash as in step 3.
13. Prepare secondary antibody Add 150  $\mu$ L per slide.
14. Incubate 30-60 minutes at RT

15. Blot excess. Wash as per step 3.
- 16.\* Prepare 250 mL of 3% # hydrogen peroxide (Sigma) in 100% methanol (Fisher). Quench endogenous peroxidases by washing for 30 min at RT. (\* optional step; use if high background level of endogenous peroxidases) (# try 0.3% H<sub>2</sub>O<sub>2</sub> first)

## **IMMUNOHISTOCHEMISTRY PROTOCOL: DLX-2 (TES-1)**

08.93

### **Protocol B: Biotin / Streptavidin / Peroxidase**

#### **Day 2 (continued)**

17. Blot excess. Wash as in step 3.
18. Prepare tertiary molecule:
19. Add 150  $\mu$ L to slide. Incubate 30-60 min at RT.
20. Blot excess. Wash as in 3.
21. Prepare DAB substrate: Waste and excess DAB treated with 2% bleach prior to discarding.
22. Add 150  $\mu$ L . Allow staining for 1- 30 minutes; usually 2-10 min. Observe under microscope.
23. Stop reaction in water for 5 min.
24. Dehydration in graded ethanols then xylene. (refer to cresyl violet protocol).
25. Mount with permount.

