

# Whole-mount in situ Hybridization

## Supplement for simultaneous detection of differently labeled riboprobes

8/15/94 Kenji Shimamura

1. Prepare digoxigenin-labeled and fluorescein-labeled riboprobes.  
(digoxigenin- and fluorescein-labeled UTP are available from Boheringer Mannheim)
2. Purify the transcripts by precipitation with lithium chloride and ethanol, or by spin-column with Sephadex G-50 beads.
3. Hybridize with the mixture of probes (0.5-1  $\mu\text{g}/\text{ml}$  each).
4. Wash probes.
5. Incubate with pre absorbed anti-fluorescein alkaline phosphatase conjugated antibodies (Boheringer-Mannheim) or anti-digoxigenin antibodies.
6. Wash antibodies.
7. Incubate with \*BCIP/NBT substrate, or Vector black alkaline phosphatase substrate (Vector Laboratories).
8. Wash with CMFET.
9. Inactivate enzyme completely by heating (70°C, 30 min), by dehydration through methanol, or by incubation with 0.1M glycine-HCl pH2.2 for 10 min.
10. Wash with TBST.
11. Incubate with the other antibodies.
12. Develop color reaction with the other set of substrate.

\*BCIP/NBT gives dark blue signal, whereas Vector black substrate gives dark brown staining.

INT/BCIP that produces brilliant red precipitates is available from Boheringer.