

Livesey lab – Oligo Array Hybridisation Protocol

Questions about this protocol should be directed to James Smith, Livesey lab, Gurdon Institute, University of Cambridge; j.smith@gurdon.cam.ac.uk

We routinely use SMART-amplified double-stranded cDNA for our array hybridisations. That cDNA is generated according to the manufacturer's instructions and used directly in the labelling reaction – we have found little improvement in labelling after cleaning up the PCR-amplified cDNA.

Hybridisation

CY DYES ARE LIGHT SENSITIVE

Labelling

Make up to 35ul DEPC H ₂ O – usually	19ul
Random primers (hex) 1ug/ul	2ul
SMART product – DNA ~400ng/ul	<u>15ul</u>
	36ul per DNA (2 DNA per array)

Denature 95C/ 5min

Snap cool on ice

1 step

Rnd Prm Bf (10x Tris, MgCl ₂ , DTT,BSA)	5ul
3NTP (0.5mM)	5ul
Klenow 10U/ul	3ul
Cy3(pink) or Cy5(purple) dye(cy3 for ref)	1ul
(1mM)	_____
	15ul

	50ul sample

Incubate 37C/ 1hour - cover with foil

2 step

Rnd Prm Bf (10x Tris, MgCl ₂ , DTT,BSA)	5ul
3NTP (0.5mM)	5ul
Klenow	1.5ul
Cy3(pink) or Cy5(purple) dye	1ul

Incubate 37C/ 30mins - cover with foil
Denature 95C/ 5min
Snap cool on ice
Add 1.5ul more Klenow
Incubate 37C/ 30mins more - cover with foil

Reaction clean-up

Shake ChromaSpin columns to re-suspend. Number columns and collecting tubes.

Break off bottom first, then lid, put into collection tube, drop in 15ml tube.

Spin at 700g for 5 mins.

Change collecting tube.

Put in the 50ul from each labelled sample to middle of gel. Spin at 700g for 5 mins.

Probes can be stored at 4C - wrap in foil.

1. Wash slides in 500ml SDS/SSC. Cover and leave for 5mins.
Remove rack, discard solution. DON'T let slides dry.
2. Cover and leave in fresh SDS/SSC again for 5mins.
Remove and blot, change rack.
3. Wash in 0.2x SSC for 1 minute
4. Remove and blot, and wash in fresh dish of 0.2x SSC for 1 minute.
5. Remove and blot, and wash in 0.1X SSC for 1 minute.

Spin at 500rpm/ 5min to dry

Rinse cover slips in 96% EtOH individually using forceps, and place in oven to dry between 2 sheets of paper towel.

Hybridisation Buffer (2ml)

1ml Formamide
0.5ml 20xSSC
0.2ml 50x Denhardt Solution
0.1 ml 10% SDS
50ul 0.2M KPhos (monobasic)
150ul DEPC
